# Sponsor Report

This section is dedicated for sponsor responsibility articles.

# Unknown Profiling of Drinking Water Using High Resolution LC-MS/MS and New Software

Caroline Ding<sup>1</sup>, Anastasia Kalli<sup>1</sup>, Charles Yang<sup>1</sup>, Tim Stratton<sup>1</sup>, Hans Grensemann<sup>2</sup>

<sup>1</sup>Thermo Fisher Scientific, San Jose, CA, USA; <sup>2</sup>Thermo Fisher Scientific, Bremen, Germany

# ABSTRACT

**Purpose:** Unknown compound profiling using high resolution LC-MS/MS and new software to confidently and quickly identify unknown compounds.

**Methods:** A treated tap water sample was collected from a city in China and stored in a plastic bottle before analysis. Mobile phase (5mM ammonium formate with 0.1% formic acid in water) was used as a blank to generate a data-dependent exclusion list for acquisition and background removal during data processing. LC-MS analysis was performed on the blank and water samples in positive modes with two replicate injections. The MS analysis employed a 70k HRMS full scan followed by top 10 data dependent ms2 collected on a Q Exactive mass spectrometer. Data analysis was performed with Thermo Scientific<sup>™</sup> Compound Discoverer<sup>™</sup> software using a single unknown data processing workflow.

**Results:** Data was processed using Compound Discoverer with one single workflow. The processing workflow included automatic unknown component detection, unknown elemental composition, library searching against mzCloud<sup>™</sup> HRAM fragmentation library, ChemSpider database search, mass list search against built-in HRAM EFS library, automatic blank removal and structure interpretation using Custom Explanations and FISh Scoring on the fly. Batch searching against mzCloud online fragmentation library proved to be the most productive and confident way for unknown compound identifications. ChemSpider database search provided more hits which complements mzCloud search, however there were too many false positives from ChemSpider search. ChemSpider search using predicted formula helped reduce the number of false positives. Built-in FISh in Custom Explanations was used to verify hits from ChemSpider against MS2 data. All in all, Compound Discoverer software provides an effective and complete workflow for unknown identifications.

# INTRODUCTION

Unknown compound profiling of water sample is very challenging due to complexity of contaminants in the water sample. Multiple software and lots of manual interpretations are usually required to identify the unknown compounds. New emerging software and tools shed light on unknown compound identifications. mzCloud is a new online HRAM fragmentation library which contains highly curated MS/MS and MSn spectra from different collision types and collision energies. It provides the fastest and most confident small molecule unknown compound ID. mzCloud search is integrated into Compound Discoverer 2.0 along with other tools like predicted compositions based on high resolution full ms, ChemSpider search that help partially identify the unknowns. This study demonstrates a simple yet powerful workflow for unknown compound profiling using high resolution Thermo Scientific<sup>™</sup> Orbitrap<sup>™</sup> mass spectrometer and Compound Discoverer software.

### MATERIALS AND METHODS

#### Sample Preparation

Untreated tap water and mobile phase blank (injection volume 50  $\mu$ L) were directly injected on to the column for chromatographic separation and MS analysis.

#### Analytical Method for LC-MS analysis

Chromatographic separation was performed with a Thermo Scientific<sup>™</sup> Dionex<sup>™</sup> Ultimate<sup>™</sup> 3000 RS LC system using a Thermo Scientific<sup>™</sup> Accucore AQ column (100 x 2.1 mm, 2.6 µ particle size). Mobile phase A was 5 mM ammonium formate and 0.1% F.A in water. Mobile phase B was 5 mM ammonium formate and 0.1% F.A in MeOH. Mass spectrometric analysis was performed on a Thermo Scientific<sup>™</sup> Q Exactive<sup>™</sup> Plus Orbitrap mass spectrometer operated in full MSddMS2 mode. Analysis was performed in positive ion mode followed by top10 data-dependent MS/MS scans. Resolution for the full MS scan was set at a 70,000 and at 17500 for the ddMS2 scans. Stepped collision energy was used at 20, 40 and 60. Analysis time, including column equilibration, was 25 min.

Retention time (min)	Flow (ml/min)	%В
0.00	0.300	2.0
0.50	0.300	2.0
2.00	0.300	40.0
20.00	0.300	95.0
22.00	0.300	95.0
22.10	0.300	2.0
25.00	0.300	2.0

# Data Analysis

Blank n/a 1	F3: Blank_pos_10ul_01
Blank n/a 2	F4: Blank_pos_10ul_02
Blank n/a 2	F4: Blank_pos_10ul_02
owater	
owater Sample Tapwa	ter 1 F15: Tapwater_pos_50ul_01

Figure 1. Sample grouping.

Samples and blanks were grouped based on user defined study factor and processed together. In this study the grouping was based on defined water type "TapWater" (see Figure 1). Sample grouping was persisted into data processing and results display.

The HRAM data was processed by Compound Discoverer software using a single processing workflow (Figure 2). The workflow employed unknown compound detection followed by online ChemSpider, mzCloud<sup>™</sup> database search and local EFS HRAM compound database search. The Mark Background Compounds node hides the background compounds in the blank files from the result table.



Figure 2. Workflow tree in Compound Discoverer software.

#### RESULTS

The results review in Compound Discoverer is broken into three parts: 1) Chromatogram view which interacts with the result table; 2) Mass Spectrum view which also interacts with the result table and displays the spectral tree for selected compound; 3) Result tables: the most important table is the Compounds table on the far left (See Figure 3). All the views can be docked, repositioned or dragged onto a second monitor.



Figure 3. Result View in Compound Discoverer.

#### **Result Filtering**

185 compounds with unique molecular weight and retention times above 1e6 peak intensity were detected from the positive mode data by Compound Discoverer excluding compounds found in the blanks. Without blank removal, the number of detected MW and RT was 711. Result filters were used to filter out compounds from the table based on user defined conditions, i.e. area threshold (see Figure 4).

#### Unknown ID with mzCloud

Compound Discoverer 2.0 includes batch compound ID against mzCloud online HRAM fragmentation library which contains high quality curated MS/MS and Msn spectra. The search algorithm allows match with ion activation energy with user adjustable ion activation energy tolerance window (Figure 5). The sophistication in the search algorithm increases the confidence in the identifications for small molecules where fragmentation pattern changes with ion activation energies.

mzCloud hits are indicated in the Compounds table with number of hits and best match scores (Figure 6). For each hit, the spectra comparison between the query spectrum and library match spectrum is visualized in a mirrored plot (Figure 7).



۵	1. Search Settings									
	Compound Classes	All								
	Match Ion Activation Type	True								
	Match Ion Activation Energy	Match with Tolerance								
	Ion Activation Energy Tolerance	20								
	Apply intensity threshold	True								
4	2. Compound Annotation									
	Assign Component Names	True								
	Assignment Threshold	90								

Figure 4. Result Filters.

Figure 5. mzCloud node settings.

	Compo	ounds 😵	Compou	inds per File	Merged Featu	ires Feature	Custom	Explanations	mzCloud Re	esults Cher	nSpider Results	Mass L	ist Search Results				
	判	Checked	Name		Pred	licted Formula		Molecular Weigh	nt RT [min]	Area (Max.)	# ChemSpide	r Results	# mzCloud Resu	ts mzCloud Best Ma -	Mass List Matches 🛨	Area	+
1	4	1	Incyclazole		C9 H	17 N3 S		189.0360	5.138	68336	0	1		2 98.4			830 650
2	4		Melamine		C3 H	16 N6		126.0653	7 1.085	247604	5	2		2 97.1			20el.48el
3	+		1-Methylben	zotriazole	C7 H	H7 N3		133.0639	133.06399 4.502		6	11		4 95.6	m		.94e1.02e1
4	-	1	Quinclorac		C10	H5 CI2 N O2		240,96999	4.565	41218	4	4		1 95.0			169el 12el
5	4		Atrazine		C81	114 CI N5		215.0942	3 7.622	21280	1	20	1	2 94.6			11e:13e
б	+		Tris(2-chloro	ethyl) phospha	te C6 F	112 CI3 O4 P		283.9542	5 6.145	82542	4	3	14	1 94.5	0		25e1.80e
7	+		PEG n5		C10	H22 O6		238.14194	4 2.952	16253	6	2		2 92.5			.63c
8	=	0	Crotamiton		C13	H17 N O		203.13110	9.739	16783	7	20		2 91.8			.68e'.59e
Pre	dicted	Checked	AMass [Da]	AMass (ppm)	Scan Number	mzCloud ID	Formula	Structure	arch Resul	ts N	lolecular Weight	Match *	Best Match *	Name			
1	0		-0.00046	-2.13	2869	42	C8 H14 CI I	CI N II	'NA		215.09377	94.6	94.6	Atrazine			
2	4	0	0.00040	1.84	2869	<u>1731</u>	C13 H13 N		NH:	Он	215.09463	21.4	21.4	3-(2-Naphthyl)-1-alanii	ne		

Figure 6. mzCloud hits in Compound Discoverer with match scores.



Figure 7. mzCloud hit spectral comparison between query and library spectra.

#### Unknown Compound Formula Prediction

Formula predictions by the Predict Composition node are listed in the sub table for each compound. The one with the best SFit% and most number of matching isotopes (#MI) is listed on the top with information like delta mass ppm. (See Figure 8)

	Co	mpo	unds 😽	Compounds per File	Merg	ged Features	Features	Custom	Explanatio	ns ma	Cloud Re	sults	ChemSp	ider Resul	ts Mass	List Sea
	F		Checked	Name		Predicted	Predicted Formula			Iolecular Weight RT		Area (Max.)		# ChemSpider F		s # mzC
1	-	=		Tricyclazole		C9 H7 N3	C9 H7 N3 S			9.03607	5.138		683360			
2	-	•		Melamine	C3 H6 N6	C3 H6 N6			5.06537	1.085		476045			1	
3	+	•		1-Methylbenzotriazole	C7 H7 N3	C7 H7 N3			3.06399	4.502		202306		i l		
4	-	- [		Quinclorac		C10 H5 C	12 N O2		240	0.96999	4.565		412184			
5	+	•		Atrazine		C8 H14 C	I N5		215	5.09423	7.622		212801		20	)
6	4	=		Tris(2-chloroethyl) phosph	ate	C6 H12 C	13 O4 P		283	3.95425	6.145		825424	1		\$
7	-	•		PEG n5		C10 H22	06		238	3.14194	2.952		162536	j		2
8	1	=		Crotamiton		C13 H17	NO	203.13			9.739 167837		167837	7		
	Hid	e Rel	lated Tab	les												
Pr	edic	ted (	Composit	ions Compounds per File	m	zCloud Results	ChemS	oider Result	ts Mass	List Sea	ch Result	s				
	F	3	Checked	d Formula		Molecular We	ght Is know	wn to Chen	Spider Δ	Mass [Da	] AMass	[ppm]	SFit [%] *	# MI •	RDBE	
3	1	-		C10 H5 CI2 N O2		240.969	973	Х		0.00026		1.07	7	5 5	8.0	
	2	-		C7 H6 Cl2 F N O3		240.970	88		-0.00089			-3.67 4		5 5	4.0	
	3	-12		C5 H8 CI2 F3 N S		240.970	240.97066			-0.00067		-2.77 4		5 4	0.0	
- 0	4	4		C6 H13 Br N P S		240.968	397			0.00102		4.24		0 4	1.0	
	5	-12		C4 H8 Br N3 O4		240.969	982			0.00017		0.72		1 3	2.0	
	6	-12		C3 H5 CI2 N7 S		240.970	042			-0.00043		-1.78		7 3	4.0	
	7	-		C5 H6 CI3 N5		240.968	240.96888		0.00111			4.62		7 3	4.0	
	8	-12		C9 H5 CI N O3 P		240.969	240.96956			0.00043		1.80		7 2	8.0	
	9	-		C7 H4 N3 O3 P S		240.97	110		-0.00111			-4.59		9 2	8.0	
	10	-		C8 H5 N O4 P2		240.969	0.96938			0.00061 2.53		2.53	1	9 2	8.0	

Figure 8. Predicted Compositions for each unknown compound listed in the sub table.



Spectral fit is visualized for each composition prediction in the spectrum window. (See Figure 9).

**Figure 9.** Spectral fit for predicted composition  $C_{10}H_5CI_2NO_2$  based on resolution.

#### ChemSpider Search and Custom Explanations

For the compounds that did not have match from mzCloud, ChemSpider hits were reviewed. ChemSpider search was performed using predicted formulas. If formula was not available, then accurate mass was used. The databases used were ACToR: Aggregated Computational Toxicology Resource; DrugBank; EAWAG Biocatalysis/Biodegradation Database; EPA DSSTox; FDA UNII – NLM. ChemSpider hits for each compound are listed in the sub table in the order of # of references (See Figure 10).

	Comp	iounds 🌾	Compo	unds per File	Merged Feat	ures Feat	tures Custom	Explanations	mzCloud R	esults	Chems	pider Results	Mass I	ist Search Results				
	2	Checked	Name		Pre	dicted Form	iula	Molecular Weig	• RT (min)	Area	(Max.)	# ChemSpide	Results	# mzCloud Results	mzCloud Best Match	Mass List Matches (*)	Area	(*)
9	0 1				CIA	4 H29 N O2		243.2201	7 12.802		115758		7	0				160
9	1 -		??		C12	2 H23 N O4		245.1631	1 7.686		1009326		5	0			.\$4e1.30	lel (lle).01e
9	2 =	E			C4	H2 F3 O5 P	S	249.9317	0.871		4834949		2	0				.84e1.83e1
	Hide F	lelated Tabl	es															
Pr	edicted	i Compositi	ons Comp	ounds per File	mzCloud F	Results C	hemSpider Resul	Mass List S	earch Resul	ts								
	₽ª!	Checked	ΔMass [Da]	∆Mass [ppm]	CSID	Formula	Molecular Wei	ght Name							Structure		#	References *
	1 🗢		0.00094	3.84	4660	C13 H19	1 245.16405	Pinacidil							to to	.×	8	1
	2 44		-0.00040	-164	<u>86846</u>	C12 H23	1 245.16270	1,1'-Iminob	1,1 -Iminobis[3-(allyloxy)-2-propanol]						4×.	yyyyy y	z	7
	3 =		-0.00040	-1.64	2289375	C12 H23	1 245.16270	2-{(2-Hydro	oxyethyl)(3-i	nethor	(ypropyl)	mino)ethyl me	thacrylat	e	X°X	xxx.J	× s	
															100	- offan	*	

Figure 10. ChemSpider hits for each compound listed in the sub table.

The problem is how do we know if any of the ChemSpider hits is the right answer? FISh Scoring in Custom Explanations in Compound Discoverer was used to verify compound ID against MS/MS data. User proposes a structure for the compound of interest, runs FISh Scoring on the fly based on the proposed structure. FISh coverage score is calculated and fragment structures are automatically annotated. Figure 11 shows an example of FISh annotations on an unknown compound based on ChemSpider proposal.

Scan this QR code on your mzCloud app to find out what this compound is! Download the mzCloud app from your App Store on iphone or androids.





Figure 11. FISh Scoring based on proposed structure for unknown ID.

#### Identified Compounds

From the 185 compounds detected in the water sample, 16 of them were identified by mzCloud automatically (see Table 1). The most dominant identified compound is Melamine. Its peak intensity is about 2e7. Others include drugs, pesticides, herbicides and etc. The ones with match score > 90 are very confident identifications based on MS/MS spectrum match and collision energy match. The other unknown compounds without mzCloud hits were much more difficult to identify. A strategy of combing ChemSpider hits, delta ppm, predicted compositions, custom explanations and FISh Scoring were used to try to identify these unknowns. However, the false positive ID rate was very high.

		Table II em	alemn semper			y mzolouc	4		
	Name	Formula	Molecular Weight	RT [min]	Delta ppm	Area (Max.)	ChemSpider Results	mzCloud Results	mzCloud Match Score
1	Melamine	C3H6N6	126.06537	1.09	0.17	2476045	2	2	97
2	Dextromethorphan	C18H25NO	271.19405	15.70	-1.62	2170464	8	1	85
3	hydroxycoumarin	C9H6O3	162.03184	14.53	-0.89	1456203	8	4	75
4	Tris(2-chloroethyl) phosphate	C6H12Cl3O4P	283.95425	6.15	-1.33	825424	3	1	95
5	Tricyclazole	C9H7N3S	189.03607	5.14	-0.03	683360	1	2	98
6	Quinclorac	C10H5Cl2NO2	240.96999	4.57	-1.07	412184	4	1	95
7	tri-phenylphophine oxide	C18H15OP	278.08637	9.66	-1.14	390088	2	3	72
8	Triethyl Phosphate	C6H15O4P	182.07087	4.93	-0.41	300639	4	2	90
9	8-Aminooctanoic acid	C8H17NO2	159.12606	1.47	0.83	286911	24	3	42
10	Atrazine	C8H14CIN5	215.09423	7.62	-2.13	212801	1	2	95
11	1-Methylbenzotriazole	C7H7N3	133.06399	4.50	0.06	202306	11	4	96
12	Crotamiton	C13H17NO	203.13118	9.74	-0.82	167837	20	2	92
13	PEG n5	C10H22O6	238.14194	2.95	-1.26	162536	1	2	93
14	phthalate	C12H14O4	222.08961	7.94	-1.79	162193	45	4	75
15	Indole-3-acetic acid	C10H9NO2	175.06345	3.47	-0.68	147127	24	2	67
16	Proline	C5H9NO2	115.06351	2.50	-1.53	105885	33	2	51

Table 1. Unknown compounds identified by mzCloud

# CONCLUSIONS

- Compound Discoverer 2.0 provides a single software solution for HRAM data processing and confident unknown compound identifications
- Unknown compound ID via batch search against mzCloud online HRAM fragmentation library proved to be the most productive and confident way for unknown compound identifications ChemSpider search combined with calculated formula from high resolution Orbitrap data complements mzCloud search but has too many false positives
- Structure elucidation using Custom Explanations and FISh Scoring in Compound Discoverer was handy and a nice way to verify ChemSpider hits against MS/MS data
- Quantitation of unknown contaminants is not the focus of this study. However, results can be exported from Compound Discoverer to software like TraceFinder for absolute quantitation.

This Sponsor Report is the reponsibility of Thermo Fisher Scientific.